


# **E-xtreme**



# **Endo Motor USER MANUAL**

Changzhou Sifary Medical Technology Co., Ltd.

Version: S02 IFU- 6035009/S02 Issued: Dec. 10 2020 Size: 160mm X 92mm


#### Content

1. Scope of E-xtreme5
1.1Parts Identification 5
1.2Components and Accessories 6
1.3 Options 6
<ol> <li>Symbols used in the User Manual7</li> <li>Before Use9</li> </ol>
3.1Intended Use 9
3.2Contraindications 9
4. Installing the E-xtreme 11
4.1Installation of the contra angle11
4.2Install the file 11
4.3Charging12
4.4Install the Insulation Sleeve13
5. Use Interface 14
5.1Panel key 14
5.2Screen display 15
5.3Terms and definition 16
6. Setting 17
6.1Set memory mode 17
6.2Set parameters······18
6.3Preset programs21
6.4Advanced setting23

6.5Parameter logic	26
7. Operation	28
7.1 Charge	
7.2Motor Operation	
8. Cleaning, Disinfection and Sterilization	
8.1Foreword	
8.2General recommendations	
8.3Autoclavable Components	
8.4Disinfection Components	
9. Error warnings	
10. Troubleshooting	40
11. Technical Data	42
12. EMC Tables	44
13. Statement	

## 1. Scope of E-xtreme

#### **1.1 Parts Identification**




- 1. Contra Angle
- 2. Motor Handpiece
- 3. Handpiece Base
- 4. Insulating Sleeve (optional)
- 5. Spray Nozzle
- Note: This product does not contain root canal file 6. Adapter

## **1.2 Components and Accessories**

Motor Handpiece	Handpiece Base	Contra Angle
(1pcs)	(1pcs)	(1pcs)
Adapter (1pcs)	Spray Nozzle (1pcs)	USER MANUAL
		(1pcs)
Certificate (1pcs)	Warranty card (1pcs)	

# 1.3 Options



# 2. Symbols used in the User Manual

	i
WARNING	If the instructions are not followed properly, operation may lead to hazards for the product or the user/patient.
ΝΟΤΕ	Additional information, explanation of operation and performance.
SN	Serial number
REF	Catalogue number
<b>**</b> *	Manufacturer
	Date of manufacture
LOT	Batch number
	Safety class II device
<b>*</b>	Type B applied part
CE	CE marking
	Direct current
	Dispose of in accordance with the WEEE directive
Ť	Keep dry
134°C	Sterilizable in a steam sterilizer (autoclave) at the temperature specified
EC REP	Authorized Representative in the European Community
-20°C -55°C	Temperature limitation
20%	Relative humidity limitation

## 2 Symbols used in the User Manual

1 06kPa	Atmospheric pressure limitation
	Manufacturer's LOGO
8	Consult instructions for use
	Be careful! Refer to relevant documents
	Washer-disinfector for thermal disinfection

#### 3. Before Use

#### 3.1 Intended Use

Use for dental root canal treatment using endodontic instruments in torque controlled continuous rotation and in reciprocating movement.

This device must only be used in hospital environments, clinics or dental offices by qualified dental personnel and not used in the oxygen-rich environment.

#### 3.2 Contraindications

The E-xtreme is contraindicated in cases where patient/user carry medical implants such as pace makers or cochlear implants etc.

Do not use the device for implants or other non-endodontic dental procedures.

Safety and effectiveness have not been established in pregnant women and children.

# WARNING

Read the following warnings before use:

1. The device must not be placed in humid surroundings or anywhere where it can come into contact with any type of liquids.

2. Do not expose the device to direct or indirect heat sources. The device must be operated and stored in a safe environment.

3. Do not use the equipment in the presence of free oxygen, anesthetic gas or combustible materials. The equipment must be operated, used and stored in a safe environment.

4. The device requires special precautions with regard to electromagnetic compatibility (EMC) and must be installed and operated in strict compliance with the EMC information. In particular, do not use the device in the vicinity of fluorescent lamps, radio transmitters, remote controls and do not use this system near the active HF Surgical Equipment in the hospital. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the E-xtreme, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

5. Please do not charge, use or store this equipment at high temperature. Please pay attention to the use and storage conditions. 6. Gloves and a rubber dam are compulsory during treatment.

7. Never open or repair the device yourself, otherwise, void the warranty.

8. If irregularities occur in the device during treatment, switch it off. Contact the agency.

9. Please use the original power adapter when charging; do not use the equipment for treatment during charging.

10. If liquid flows out of the handpiece, it can be judged as battery leakage. Please stop using immediately and contact the local dealer for treatment.

11. Do not dismount the contra angle during the operation of the main engine, otherwise the contra angle and motor gear will be damaged.

12. Please use the original contra angle.

13. Use continuous file in continuous mode; use reciprocating file in reciprocating mode, and use according to rotation speed, torque and return angle recommended by the root canal file manufacturer.

## 4. Installing the E-xtreme


#### 4.1 Installation of the contra angle

Make sure 4 pins on contra angle alignment the slots of handpiece, plug them together until it "clic- k" securely into place.
The contra angle can be 360 deg- rees rotated without take off, make it easy to watch the LCD in trea- tment by rotating the contra angle.

# 

 Make sure the assembly is connected properly, otherwise might cause unexpected motor reverse, even hurt the patients
 After connecting the contra angle and handpiece, pull it gently to make sure the connection is good.

## 4.2 Install the file



## 

1. Inspect the file head before inserting the file. Do not use the damaged file head.


2. Be careful when inserting and removing files to avoid injury to fingers.

3. Pull the file gently to make sure that the file is secure in handpiece properly, otherwise it may pop out and hurt the patient.

4. When removing the file, press the button tightly to release the internal bayonet. If the bayonet is not fully released, the bearing will be damaged.

5. Make sure the motor is stopped when inserting and removing files.

# 4.3 Charging





# 4.4 Install the Insulation Sleeve

R	Install: assemble according
	to the left figure
	Remove: pull out the cover
	in the opposite direction
	The insulating sleeve is ma-
	inly used for secondary isol-
	ation to avoid cross infec-
	tion.

#### 5. Use Interface

#### 5.1 Panel key






#### 5.2 Screen display

## 5.3 Terms and definition

Fwd/Fw	Forward ( Clockwise rotation )
Rev/Rv	Reverse ( Counter clockwise rotation )
REC	Reciprocation: Be applied to reciprocating file, path file and rotary file protection by setting some special angle
Memory mode	Such as M0-M9
Operation mode	Such as Fwd, Rev (set in M1-M9), Reciprocation (M0)
TRQ	Torque
MEM	Memory
R·D	Rotate Direction
DIR	Direction
Separation of instruments	The file used in root canal therapy is broken accidentally.

### 6. Setting

#### 6.1 Set memory mode




### 6.2 Set parameters

# WARNING

All parameters must be set according to the recommended values of root canal file manufacturer. Before starting the device for operation, make sure that all parameters are correct, otherwise there is a risk of instrument separation.

	Before starting the motor, check
	whether the operation mode (1)
	is correct. If it is not the expected
	operation mode, press S once
M 1 3.0 N·cm	on the standby mode to enter
	the operation mode selection,
	press + or - to switch, and
	then press • or wait for 5 se-
	conds to confirm the operation
	mode.
••••• 400 RPM <b>( )</b>	The left figure describes the set-
M 1 3.0 N·cm	ting of common functions of the
	device (memory mode M1 to M9
Jer	are applicable). In standby mode
• SPEED: 4002	(1)press S to enter the speed
	setting, press + or - to set
S	the speed value; after the speed
	value is set, press S again to
$\begin{array}{c} \text{IIII}  \text{TORQUE: } 3.0 \\ \text{M 1}  \text{R-D}  \text{IIII}  \text{IIIII} \end{array}$	enter the torque setting, press
	+ or - to set the torque va-
Set.	lue. After the torque value is set,
IIII DIR: Fwd	press S again to enter the ro-
	tation direction operation mode
SE	(two modes: Fwd and Rev), press
	+ or - to set, then press •
400 RPM	or S or wait for 5 seconds to
	confirm the operation mode.

Page 18 / 52




6 Setting

1 2 DIR: Fwd M 1 MEM FW RV 3 4	Fwd or Rev can be set in opera- tion mode of M1-M9. Press $S$ till the operation mode (124) is displayed, press + or - to set, then press • or wait 5 seconds to confirm. The left figure(3) means press $S$ ag- ain will enter memory mode (MEM). Mem Rev operation mode is se- lected, a continuous alarm sound will appear after the device is started to remind the operator that the motor is in reverse rotation.
$\begin{array}{c} 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 \\ 1 $	The left figure describes the set- ting of M0(REC mode):Press <b>S</b> at M0(①)to switch different reci- procating preset programs. Th- ere are altogether 5 sets of pre- set reciprocating programs.

## 6.3 Preset programs

1 M1 Fwd Protaper SX&S1 400RPM Protaper S2 3.0N • cm 2	For the convenience of the operator, so- me common root canal file systems are preset. Long press S at standby mode to en- ter preset mode, the screen will display as the figure shows on the left. M1(①)means current memory mode, o- perator can choose preset mode(②) to replace it. Press + or - to switch (press the key several times to jump to the next page). Then press o to con- firm.
One Shape One Flare One Curve 1 4 3 3	If you select one of the preset program- ms, such as "one curve" (①), the corr- esponding operation mode(②), speed (③) and torque(④) will be automatically set. NoTE Protaper <sup>®</sup> , GATES <sup>®</sup> , Pro.Glider <sup>®</sup> , and Wave one <sup>®</sup> is a registered trademark of Dentsply. Mtwo <sup>®</sup> , Flex.Master <sup>®</sup> , Reciproc <sup>®</sup> and R- Pilot <sup>®</sup> is a registered trademark of VDW. K3XF <sup>®</sup> , TF <sup>®</sup> is a registered trademark of SybronEndo. OneG <sup>®</sup> , OneShape, OneFlare, 2Shape and OneCurve <sup>®</sup> is a registered tradem- ark of Micro-Mega XPendo.Shaper <sup>®</sup> , XPendo.Finisher <sup>®</sup> , Ir- ace <sup>®</sup> , BT-Race <sup>®</sup> and BioRace <sup>®</sup> is a reg- istered trademark of FKG



When the preset mode is selected, the memory number(1) will be changed to the preset name, operation mode (2), speed (3) and torque (4) will also be automatically set.

# ΝΟΤΕ

All memory modes (from M1 to M9) can be replaced by preset programs in this way.

# 6.4 Advanced setting

× SET	Hold <b>S</b> then press <b>•</b> (for about 0.5 seconds) at off mode to enter advanced set mode.
$ \begin{array}{c} 1 \\ 2 \\ \hline BEEP VOL: Low \\ \hline O \\ \hline 3 \end{array} $	It will be about 1 second at SET logo an- d then enter ① "BEEP VOL" (beep vol- ume set). Press + or - to set (②) (Mute、Low、Mid、High), then pres- s to confirm and shut down. The left figure③means press S once in t- his case will enter auto power off time set (A.P).
AUTO P.W.R: 10min	Press S again at beep volume set 1 will enter " AUTO P.W.R" (auto power off time set), press + or - to set 2 (3-15 minutes), then press • to confirm and shut down. The left figure 3 means press S again will enter a- uto set confirm time (S.T).
1 2 SET TIME: 5 Sec 1 2 SET TIME: 5 Sec 1 2 3	Press <b>S</b> again at auto power off time set ① will enter " SET TIME " (auto set confirm time). Press + or - to set ② (3-15 seconds) , then press • to c- onfirm and shut down. The left figure ③ means press <b>S</b> again will enter hand habit set mode (L.R).
$ \begin{array}{c} 1 \\ 2 \\ Habit Hand: Right \\ \hline 0 \\ \hline 3 \\ 4 \end{array} $	Press S again at auto set confirm ti- me 1 will enter" Habit Hand" (hand habit set). Press + or - to set24 (Left, Right), then press • to co- nfirm and shut down. The left figure 3 means press S again will enter auto calculation(CAL). After switching to the left-handed habit, the display interface will rotate 180 ° to facilitate the left-handed operator to observe the display screen.

6 Setting

r	o setting
Calibration Calibration Calibration Calibration Calibration Calibration Calibration Calibration Calibration Calibration Calibration Calibration Calibration Calibration	Press <b>S</b> again at hand habit set will enter <b>1</b> "Calibration" auto calibration, press <b>+</b> or <b>-</b> to set <b>3</b> (ON, OFF), Choose "ON" and press <b>o</b> or <b>S</b> to confirm and the device will auto cali- brate.The left figure <b>2</b> means when choose "OFF" and press <b>S</b> again w- ill enter restore settings (R.S).
Calibration -1	When the device auto calibrate①,the screen will display the progress②(thro-ugh the bar), after calibrating③ the bar will be full and the device will beep out.
	Before automatic calibration, make sure that the original contra angle is insta- lled on the handpiece, and the root canal file is not installed on the contra angle. If the contra angle is not insta- lled or non-original contra angle is installed, the torque after calibration may be incorrect, which may bring the risk of instrument separation.
Restore Setting 2 3	Press <b>S</b> again at auto calibration with swtich" off" will enter (1) "Restore Set- ting (restore settings). Press $+$ or $-$ to set(3)(ON, OFF). Choose "ON" and press • or <b>S</b> to confirm and restore all settings. Ch- oose "OFF" and press <b>S</b> again (2) (MEM) to confirm, save all settings ab- ove and return to the memory mode. All parameters will be covered by the factory default parameters. (refer to chapter 7.5 parameter logic). After restoring the factory settings, the parameters set by the customer will be covered by the factory default

6 Setting

5
parameters. If necessary, please reco-
rd the parameters before restoring th-
e factory settings.

#### 6.5 Parameter logic

The factory default parameters of the ten memory modes are shown in the table below. The parameters can be adjusted as needed.

The default advanced settings parameters are shown in the following table. The parameters can be adjusted as needed


Parameter	M0	M1	M2	M3	M4	M5	M6	M7	M8	M9
Operation mode	REC	Fwd	Fwd	Fwd	Fwd	Fwd	Rev	Rev	Fwd	Fwd
Speed(rpm)	N/A	350	300	400	400	300	350	500	500	650
Torque(N•cm)	N/A	2.5	3.0	2.0	1.5	1.5	2.5	2.0	2.5	2.0
Fwd angle	30	N/A								
Rev angle	150	N/A								

The default advanced settings parameters are shown in the following table. The parameters can be adjusted as needed.

Volume BEEP VOL	Mid
Auto power off AUTO P.W.R	10min
Auto set confirm SET TIME	5s
<b>Hand habit</b> Habit Hand	Right

Auto calibrate Calibration	Off
Restore settings Restore Setting	Off
/	/
/	/

**The rotational speed (RPM)** varies in different operating modes, as shown in the table below



**Torque** (N•cm) in different operation modes, the torque value can be set differently even in the same operation mode when the speed value is set differently. See the table below for details.

	•	
Fwd	Rev	REC
0.50.81.01.53.23.54.0	1.8 2.0 2.2 2.5 3.0	/

There are 5 fixed values of **reciprocating Angle** in M0 reciprocating mode, and the Angle cannot be changed., as shown in the table below.


$\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{$	Fwd	Rev	REC
			Five sets of fixed values
			1. Fwd angle 30°, Rev angle 150°
reciprocating	/		2. Fwd angle 150°, Rev angle 30°
Angle		3. Fwd angle 180°, Rev angle 30°	
			4. Fwd angle 210°, Rev angle 30°
			5. Fwd angle 250°, Rev angle 30°

# 7.Operation

# 7.1 Charge


	Display the remaining power. The remaining charge is less than 15%.
	<ul> <li>NOTE</li> <li>1. If the battery power is less than 15%, it must be recharged within 30 days, otherwise the battery will be irretrievably damaged due to low power.</li> <li>2. If you do not use this product for a long time, please charge the product at least once a month.</li> </ul>
<b>Low Power</b> Please Charge	If the battery power is lower than 15%, the speed and to- rque may be lower than the set value. As shown in the left fig- ure, the low power alarm will appear on the display screen with continuous use, and the device will automatically shut down.
	<b>D</b> Because the display of the rema- inning power is based on the voltage level, if a large torque load appears during the oper- ation, the remaining power display may appear a short-term decrease.

7 Operation



#### 7 Operation

#### 7.2 Motor Operation



In the standby mode, the root canal preparation device is started by pressing the main switch , After startup, the progress bar will be displayed on the display screen (for details of the progress bar, please refer to chapter 5.2 display screen interface).



1. Before using in the treatment, please try it out of the mouth to ensure that the function of the device is normal.

2. The root canal file may be damaged suddenly when it enters into the root canal which is too curved or not in good shape. When the user feels that the root canal is abnormal, please stop using the device immediately and confirm the correct operation parameters and methods.

3. Even if the normal parameters are set, due to the metal fatigue of the root canal file, the instrument will be separated. Therefore, when using the root canal file, do not exceed the times recommended by the manufacturer, and replace it in time.

4. When the root canal file is subjected to excessive external force, it may break. When using this equipment, do not apply excessive external force to the root canal file.

5. Do not press the back cover of the contra angle during the treatment, otherwise the equipment will be damaged, and even the flying file will hurt the patient.

6. The electromagnetic noise in the surrounding environment may interfere with the normal operation of the equipment. Please do not completely rely on the automatic control of the equipment, and always pay attention to the feedback information on the LCD screen.

# 

1. When there is any abnormality, please stop using the equipment. This equipment is not suitable for all types of root canals. It is recommended to use according to the instructions of root canal file.

2. The root canal file is easy to fracture at high speed. Please follow the rotation speed recommended by the manufacturer. Please check the set speed before use.

3. When using this equipment, use the root canal file with materials other than nickel titanium carefully.

4. Please use disposable gloves and rubber barrier for treatment.

5. After the treatment, please take out the root canal file to avoid damage to the root canal file.

#### 8. Cleaning, Disinfection and Sterilization

#### 8.1 Foreword

For hygiene and sanitary safety purpose, the components (contra angle, and insulating sleeve) must be cleaned, disinfected and sterilized before each usage to prevent any contamination. This concerns the first use as well as the subsequent uses. Comply with your national guidelines, standards and requirements for cleaning, disinfection and sterilization. Reprocessing procedures have only limited implications to this dental device. The limitation of the numbers of reprocessing procedures is therefore determined by the function / wear of the device. From the processing side there is no maximum number of allowable reprocessing. The device should no longer be reused in case of signs of material degradation.


In case of damage the device should be reprocessed before sending back to the manufacturer for repair.

#### 8.2 General recommendations

• The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments, where applicable after sterility.

- For your own safety, please wear personal protective equipment (gloves, safety glasses, etc.).
- Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, and FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer.
- The water quality has to be convenient to the local regulations especially for the last rinsing step or with a washer-disinfector.
- Thoroughly clean and wash the components before autoclaving.
- Do not lubricate the motor handpiece.
- Do not clean the contra angle with an ultrasonic cleaning device.
- Do not use bleach or chloride disinfectant materials.

#### 8.3 Autoclavable Components



	sing area to avoid any damage and contaminate on to the environment.
Preparation for Decontaminati- on	<ul> <li>The devices must be reprocessed in a disassembled state.</li> <li><i>WARNING</i></li> <li>1. Do not fail to take out the file before cleaning the contra angle.</li> <li>2. Observe suitable personal protective measurees.</li> </ul>
Pre-Cleaning	Do a manual pre-cleaning, until the components are visually clean. Submerge the components in a cleaning solution and flush the lumens with a water jet pistol with cold tap water for at least 10 seconds. Clean the surfaces with a soft bristol brush.
Cleaning	Regarding cleaning/disinfection, rinsing and drying, it is to distinguish between manual and automated reprocessing methods. Preference is to be given to automated reprocessing methods, especially due to the better standardizing potential and industrial safety. Automated Cleaning: Carefully put the components into the washer-disinfector on a tray and set the parameters as follows, then start the program: • 4 min pre-washing with cold water (<40°C); • emptying • 5 min washing with a mild alkaline cleaner at 55°C; • emptying • 3 min neutralising with warm water (>40°C); • emptying


	• 5 min intermediate rinsing with warm water
	(>40°C);
	• emptying
	<i>The automated cleaning processes have been validated by using 0.5% neodisher MediClean forte (Dr. Weigert).</i>
	Note Acc. to EN ISO 17664 no manual reprocessing methods are required for these devices. If a manual reprocessing method has to be used, please validate it prior to use.
	<ol> <li>Use only approved washer-disinfectors according to EN ISO 15883, maintain and calibrate it regularly.</li> </ol>
	<ol> <li>Follow instructions and observe concentrations given by the manufacturer (see general recommendations).</li> </ol>
	3. Avoid any contact between the contra angle and any instrument, kit, support or container.
	Automated Thermal Disinfection in washer/disinfector under consideration of national requirements in regards to A0 value (see EN ISO 15883).
Disinfection	A disinfection cycle of 5 min disinfection at 93°C has been validated for the device to achieve an A0 value of 3000.
	After manual cleaning, the instruments should be automated disinfected of sterilized immediately. A manual disinfection is not recommended.
Drying	Automated Drying: Drying of outside of instrument through drying cycle of washer/disinfector. If needed, additional manual drying can be performed through lint free towel. Insufflate cavities of instruments by

	using sterile compressed air.
Functional	Visual inspection for cleanliness of the
Testing,	instruments and reassembling. Functional testing
Maintenance	according to instructions of use. If necessary, perform reprocessing process again until instrument is visibly clean. Before packaging and autoclaving, make sure that the components have been maintained acc. to manufacturer' s instruction. Only the contra angle needs to be lubricated.
	Black oil
	1. Before autoclaving, the contra angle must be lubricated.
	<ol> <li>Attaching the spray nozzle to oil can and contra angle, press the oil can button more than 3 seconds, till all the black oil flow out from the head of the contra angle.</li> </ol>
Packaging	Pack the instruments in an appropriate
	packaging material for sterilization.
	<ol> <li>Check the validity period of pouch given by the manufacturer to determine the shelf life.</li> </ol>
	<ol> <li>Use pouches which resist to a temperature up to 141°C and in accordance with EN ISO 11607.</li> </ol>
Sterilization	Sterilization of instruments by applying a fractionated pre-vacuum steam sterilization process (according to EN 285/EN 13060/EN ISO 17665) under consideration of the respective country requirements.

	Minimum requirements: 3 min at 134 °C (in EU: 5 min at 134 °C)
	Maximum sterilization temperature: 137°C Flash sterilization is not allowed on lumen instruments!
	1. Use only approved autoclave devices according to EN 13060 or EN 285.
	2. Use a validated sterilization procedure according to EN ISO 17665.
	3. Respect the maintenance procedure of the autoclave device given by the manufacturer.
	<ol> <li>Use only this recommended sterilization procedure.</li> </ol>
	5. Control the efficiency (packaging integrity, no humidity, color change of sterilization indicators, physicochemical integrators, digital records of cycles parameters).
	6. The sterilization procedure must comply with EN ISO 17665.
	7. Wait for cooling before touching.
Storage	Storage of sterilized instruments in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.
	1. Sterility cannot be guaranteed if packaging is open, damaged or wet.
	2. Check the packaging and the contra angle before using it (packaging integrity, no humidity and validity period).

The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for use. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

### 8.4 Disinfection Components



## 9.Error warnings

<b>Overload</b> Restart Motor	This warning will appear on the display screen if the load exceeds the capacity of the standby machine during reversal. Please press the main switch key to restart the standby machine.	
<b>Low Power</b> Please Charge	The power is very low, charge it immediately.	

### 10. Troubleshooting

When trouble is found, check the following points before contacting your distributor. If none of these are applicable or the trouble is not remedied even after action has been taken, the product may have failed. Contact your distributor.

Problem	Cause	Solution	
Cannot power on Charging indication does not	power onpressing the main switch is too short.for more secondsChargingWrong adapter is usedPlease use nal power adapter		
appear on the handpie-	is not plugged into the socket The socket is not	Please check the con- nection Please check the con-	
ce screen	energized	nection	
Handpie- ce screen does not display any informati- on	The handpiece is damaged	Long press the main switch to start the device, check whether the sound is normal, and press the main switch again to check whether there is the sound of motor rota- tion. Then contact the dealer.	
The motor does not rotate	Contra angle stuck	Pull out the contra angle and check whether the motor rotates. If it can rotate normally, please clean or repair the contra angle	

	The handpiece is protected or damaged by the system	Check according to error warning	
The motor cannot stop	There is a short circuit in the internal circuit	Press the "s" key to stop the motor and contact the dealer	
The motor reverses	The reverse value of the torque setting is reached	Check whether the torque limit is too small	
uncontro- llably	Rev inversion mode is set	Check settings	
Motor does not reverse	Excessive torque reversal value is set	Check settings	
Frequent switching between forward and reverse rotation of motor	Rec (reciprocating) operation mode is set	If not, switch the opera- tion mode	
No sound	Volume set to "mu- te"	Set the volume to low, mid, or high	
The handpie- ce sends out a continuo- us alarm.	Rev is set	If the setting is expected, ignore the alarm	

## 11.Technical Data

Manufacturer	Changzhou Sifary Medical Technology Co., Ltd
Model	E-xtreme
Dimensions	17.5cm x 10.9cm x 8.4cm±1cm (Outer box)
Weight	0.6kg±15%
Contra angle	Compatible with rotary and reciprocating instruments, equipped with 2.35mm nickel titanium root canal file conforming to ISO 1797:2017, Type 1, Files length 11-31mm.
Power supply	Lithium ion battery: 3.7V, 800mAh, ±10%
Charger power supply	AC 100-240 V, ±10%
Charger power output	5V 1A
Frequency	50/60Hz, ±1Hz
Charger nominal power input	0.4A Max
Torque range	0.5 – 4.0N·cm
Speed range	120-650 rpm
Type of protection against electrical shock	Class II and internally powered equipment
Applied part	В
Operation mode	Intermittent operation, working for 60 minutes / stopping for 5 minutes
Ingress Protection	IPX0
Operation conditions	Use: in enclosed spaces Ambient temperature: 5°C ~ 40 °C Relative humidity: <80% Operating altitude < 3000m above sea level

Transport and storage conditions	Ambient temperature: -20 °C ~ +55 °C Relative humidity: 20% ~ 80 % Atmospheric pressure: 70kPa ~ 106kPa
----------------------------------	---

### **12.EMC** Tables

# Guidance and manufacturer's declaration - electromagnetic emissions

The **E-xtreme** is intended for use in the electromagnetic environment specified below. The customer or the user of the **E-xtreme** should assure that it is used in such an environment.

Emissions test	Complian- ce	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The <b>E-xtreme</b> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The <b>E-xtreme</b> is suitable for use in all establishments, incl-
Harmonic emissions IEC61000-3-2	Class A	uding domestic establishen- ts and those directly connect-
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	ed to the public low voltage power supply network that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration – electromagnetic immunity

The **E-xtreme** is intended for use in the electromagnetic environment specified below. The customer or the user of the **E-xtreme** should assure that it is used in such an environment.

lmmunity	IEC 60601	Compliance	Electromagne-
test	test level	level	tic
			environment - guidance

Electrostatic discharge (ESD) IEC	+/- 8 kV contact	+/- 8 kV contact	Floors should be wood, concrete or
61000-4-2	+/- 2 kV, +/- 4	+/- 2 kV, +/- 4	ceramic tile. If
	kV, +/- 8 kV,	kV, +/- 8 kV,	floors are
	+/- 15 kV air	+/- 15 kV air	covered with
			synthetic
			material, the
			relative
			humidity
			should be at
			least 30 %.
Electrical fast	±2kV	±2kV	Mains power
transients/b-	100kHz	100kHz	quality should
ursts IEC	repetition	repetition	be that of a
61000-4-4	frequency	frequency	typical
			commercial or
			hospital
			environment.
Surge IEC	Line to line:	Line to line:	Mains power
61000-4-5	±0.5kV, ±1kV	±0.5kV, ±1kV	quality should
			be that of a
	Line to earth:	Line to earth:	typical
	±0.5kV, ±1kV,	±0.5kV, ±1kV,	commercial or
	±2kV	±2kV	hospital
			environment.

11 Technical Data

			Mains power		
			quality should		
Voltage dips			be that of a		
IEC	0% UT; 0.5	0% UT; 0.5	typical		
61000-4-11	cycle	cycle	commercial or		
	at 0°, 45°, 90°,	at 0°, 45°, 90°,	hospital		
	135°, 180°,	135°, 180°,	environment. If		
	225°, 270°,	225°, 270°,	the user of		
	and 315°	and 315°	devices require		
			continued		
	0% UT; 1 cycle	0% UT; 1 cycle	operation		
	and 70% UT;	and 70% UT;	during power		
	25/30 cycles	25/30 cycles	mains		
	sine phase at	sine phase at	interruptions, it		
	0°	0°	is recommend-		
			ded that devic-		
Voltage	0% UT;	0% UT;	es be powered		
interruptions	250/300 cycle	250/300 cycle	form an		
IEC			uninterruptible		
61000-4-11			power supply		
			or a battery		
Rated Power	30 A/m	30 A/m	Power freque-		
frequency	50Hz or 60Hz	50Hz or 60Hz	ncy magnetic		
magnetic			field should be		
field IEC			at levels chara-		
61000-4-8			cteristic of a		
			typical location		
			in a typical		
			commercial or		
			hospital		
			environment.		
Note: UT: rated voltage(s); E.g. 25/30 cycles means 25 cycles at 50Hz					
or 30 cycles at	60Hz				

# Guidance and manufacturer's declaration – electromagnetic immunity

The **E-xtreme** is intended for use in the electromagnetic environment specified below. The customer or the user of the **E-xtreme** should assure that it is used in such an environment.

Immunity test Conducted disturbances induced by RF fields IEC 61000-4-6	IEC 60601 test level 3 V 0.15 MHz – 80 MHz, 6 V in ISM bands between 0.15 MHz and 80 MHz, 80 % AM at 1 kHz	Compli- ance level 3 V	Electromagnetic environment - guidance Portable and mobile RF communications equipment should be used no closer to any part of the E-xtreme, including cables, than the recommended
Radiated RF EM fields IEC 61000-4-3	3 V/m, 80 MHz – 2,7 GHz, 80 % AM at 1 kHz	3V/m	separation distance calculated from the equation applicable to the frequency of the transmitter.
Proximity fields from RF wireless communicate- on equipment IEC 61000-4-3	See the RF wireless communicate- on equipment table in "Recommend- ed minimum separation distances"	Compli- es	minimum separ- ation distances See the RF wireless communication equipment table in "Recommended minimum separa- tion distances"

#### **Recommended minimum separation distances**

Nowadays, many RF wireless equipments have being used in various healthcare locations where medical equipment and/or systems are used. When they are used in close proximity to medical equipment and/or systems, the medical equipment and/or systems' basic safety and essential performance may be affected. The **E-xtreme** has been tested with the immunity test level in the below table and meet the related requirements of IEC 60601-1-2:2014. The customer and/or user should help keep a minimum distance between RF wireless communications equipments and the **E-xtreme** as recommended below.

Test frequ- ency (MHz)	Band (MH- z)	Service	Modu- lation	Maxim- um power (W)	Dist- ance (m)	lmmu- nity test level (V/m)
385	380- 390	TETRA 400	Pulse modu- lation 18Hz	1.8	0.3	27
450	430- 470	GMRS 460 FRS 460	FM ± 5 kHz deviat- ion 1 kHz sine	2	0.3	28
710 745 780	704- 787	LTE Band 13, 17	Pulse modu- lation 217Hz	0.2	0.3	9
810	800-	GSM 800/900 , TETRA	Pulse modu-	2	0.3	28
870	960	800, iDEN 820,	lation 18Hz		0.0	

930		CDMA 850, LTE Band 5				
1720		GSM 1800; CDMA 1900;				
1845	1700 -199 0	GSM 1900; DECT; LTE	Pulse modu- lation 217Hz	2	0.3	28
1970		Band 1, 3, 4, 25; UMTS				
2450	2400 -257 0	Bluetoo- th, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modu- lation 217Hz	2	0.3	28
5240	5100 -580 0	WLAN 802.11 a/n	Pulse		0.3	9
5500			modu- lation	0.2		
5785			217Hz			

# WARNING

1. Use of accessories and cables other than those specified or provided by the manufacturer of **E-xtreme** could result in increased electromagnetic emissions or decreased electromagnetic immunity of **E-xtreme** and result in improper operation.

#### **Cable information:**

Cable Name	Cable Length	Shielded or	Remark	
	(m)	not		
Adapter Cable	1.2	No	/	

2. Use of **E-xtreme** adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, **E-xtreme** and the other equipment should be observed to verify that they are operating normally.