

Wide Spectrum Curing Expert

#### **GUILIN WOODPECKER MEDICAL INSTRUMENT CO., LTD.**

Registration Add & Manufactory Add: Information Industrial Park, Guilin National High-Tech Information Park 541004 Guilin, Guangxi China

Europe Sales Dept.: +86-773-5873196, +86-773-2125222

North America, South America &

Oceania Sales Dept.: +86-773-5873198, +86-773-2125123 Asia & Africa Sales Dept.: +86-773-5855350, +86-773-2125896


Fax: +86-773-5822450

ZMN/WI-10-071 V1.1

E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com Website: http://www.glwoodpecker.com

The pictures are only for reference.

#### Distributor:







#### The **New Generation** of Curing light

### **B-CURE**


ONE CURE<sup>™</sup>



# New generation of wide-spectrum LED



Three built-in LEDs in one light, which can achieve an optimal broadband spectrum of 385 to 515 nm that is similar to the spectrum of halogen lights; therefore, new light is suitable for polymerizing the entire range of photoinitiators and materials currently used in dentistry.



The light-absorbing peak of traditional material (including camphorquinone; the red line on the figure) is 460nm, so the curing can be done by either general Curing Light or Wide Spectrum Curing Light.

But camphorquinone is of a yellow color.

It is not suitable to add large amounts of it to bleaching resin or light-colored resins.

The light-absorbing peak of new materials (including TPO; the green line in the figure above) is 405nm.

Compared with general LED Curing Light,

B-Cure owns additional 385nm-420nm wave band which is able to better curing this kind of new material.

TPO is of white color. Thus it is suitable for adding to the bleaching resin and light-colored resin.

#### 3 seconds for curing

B-Cure is equipped with 5w LED.


The maximum light intensity can reach 2000mw/cm² with such a small size.

Curing normal resin that is thicker than 2mm in three seconds.

Play an important role in cementation of all-ceramic restoration and orthodontic bracket.

#### **Ortho mode**

In the ortho mode that is unique to B-Cure, the light intensity automatically adjusts to 2000mw/cm². The user can set the single curing time for 3s or 5s, and the device will automatically conduct 10 circles of curing with interval of 1s. One operation can basically finish the cementation of 2 to 3 brackets, which reduces operations and improves efficiency.




10 cycles with an interval of 1s after each cycle; Time setting for each cycle can be 3s or 5s.




#### Under Soft mode.

Under Soft mode, the power of LED gradually increases, giving a longer gel period for resin to reduce the edge discomfort caused by stress contraction of resin.



Constant high intensity of 1800-2000 mw/cm<sup>2</sup>; Time setting: 3s and 5s.



Light intensity increases from 0 mw/cm<sup>2</sup> to 1200 mw/cm<sup>2</sup> in 5 sec. Time setting: 5s, 10s, 15s and 20s.

#### Pulse mode


**Turbo mode** 

# 1200 light intensity(mw/cm²) 0 0.8 1.0 1.8 2.0 Time(s)

Continuous cycle of 1s, including the working time of 0.8s and an interval of 0.2s; Time setting: 5s, 10s, 15s and 20s.

#### **Normal mode**

Soft mode




Constant intensity of 1000-1200 mw/cm²; Time setting: 5s, 10s, 15s and 20s.

#### **Designed for comfort**

#### **Optimized Optical Fiber**

Can be autoclaved under 134°C high temperature and 0.22MPa high pressure, effectively avoiding cross infection.





#### **Continuous power supply**

The battery is easy to change by pressing the button on the device.

And the battery can be charged separately by putting it on the charger.

Two original batteries can keep the device working continuously.



#### More humanfriendly design

B-Cure adopts
ergonomic design with
smooth lines and small size.
Doctors can easily hold it.





# LED display screen, better visibility

Clear vision at a glance.


#### **Parameters**

Power input: AC100V-240V 50Hz/60Hz Light intensity:1000mW/cm²-2000mW/cm²

## Constant light intensity

Constant light output.

Solidification effect will not be affected by the consumption of remaining power.



▲ Tests have shown: curing light without constant light intensity will cause incomplete solidification to the bottom of resin when power drop to about 60%.

#### **Light hood**

Studies prove that Blue-ray will cause macular degeneration to retina and eyes.

Woodpecker light hood is made of selected quality materials, rejecting Blue-ray injury.



