

Eco 450

Features

- Automatic safety overload protection system
- Less noise, less vibration
- Can be laid on the side
- Variable speed control pedal (FS-60N)



▲ Laid on the side for narrow space



Specifications control box only

- Power supply : AC100 - 120V or 200 - 240V 50/60 Hz
- Electric power : 50W
- Speed : MAX 40,000 RPM
- Weight : 1.6 kg
- Dimensions : W141 x D186 x H97(mm)
- Contents : Control unit, Handpiece, Foot control, Handpiece stand