




LUCKCOME


L11T 2.8inch TFT Bilirubin Meter



CARE LIFE FROM THE BEGINNING



2.7inch Color touch screen




Built-in color calibration sensor




Wireless/Wired charging function



WiFi/Bluetooth data transmission



Patient ID scan



Buzzer sound during power On/off & after measurement

Features

- Color touch screen, easy to use
- Built-in color calibration sensor , supporting value 0 and 20 calibration
- Wireless charging
- Wi-Fi/bluetooth networking function
- Support code scanning
- Supports multi-user continuous measurement
- Long lifetime LED light source
- Sound indication during power on/off and completion of measurement
- Low voltage indication
- Support date and time modification
- Supports touch screen key sounds on/off
- 5-level screen brightness
- Average value is based on 1-5 times of measurement
- Measurement units: mg/dl or $\mu\text{mol/l}$
- Auto off after 1 to 5 minutes of non-operation.
- History data: Nurse No. ID No. ,Neonate ID No., measurement results & time, etc.
- 200 sets history data
- Rechargeable Li-Ion battery

Device Information

