

HDMI HD Intraoral camera

Model: MD-3000



Patented product



16:9



HDMI plug

Product features:

- HDMI HD 1080P 16:9 full display.
- Quad screen and full screen display.
- Save pictures into U disk, easy to copy pictures to computer.
- Wifi transfer pictures to computer (Support windows system only).
- Can page up/down and delete picture on the handpiece.



12 White LEDs

Specifications

Sensor	1/2.7" Sony FHD CMOS
TV standard	PAL
Resolution	1080P; 720P
Light	12 white LEDs
Focus distance	3~50mm
Angle of view	105°
Cable length	2M
Output port	HDMI (16:9)& VGA

Plugs



- * IR port
- * HDMI output port
- * VGA output port
- * U disk port
- * Power port (DC 12V)