

## PRODUCT SPECIFICATION

Product Name:	POWDER JET HANDPIECE
Model:	ProJet 2

Preview:





**Function & Introduction:**

- Function:
- Remove plaque for placement of sealants.
- Surface Preparation prior for bonding/cementation of inlays onlays, crowns and veneers.
- Surface preparation prior to placing composite restorations.
- Cleaning prior to bonding orthodontic brackets.
- Effectively remove plaque and stain for orthodontic patients.
- Cleaning implant fixture.
- Remove stain for shade determination.
- Remove plaque prior to fluoride treatment.
- Polishing on the surface of teeth after scaling treatment by piezo scaler.
- Introduction:
- Autoclavable detachable handpiece, PJ-01 supragingival sandblasting handpiece and PJ-02 periodontal sandblasting handpiece fitting with ProJet 2 with 3 levels powder adjustments.
- The Connection Type is KAVO Multiflex/ KAVO Multiflex LED.
- It suits any voldent dental powder.

**Specifications:**

Classification 93 / 42 EEC	Class II a
Water input pressure	0.5Bar~2.0Bar, (50~200kPa)
Air inlet pressure	3.5Bar~4.5Bar, (350~450kPa)
Air Flow Rate	2-15l/min
Water Flow Rate	20~100ml/min
Box Size*1 set	26.0cm*19.0cm*5.5cm
Packing Size*12 sets	41.0cm*41.0cm*28.0cm
Weight of main unit	145.0g ( ± 5% )
Gross Weight*1 set	520.0g ( ± 5% )
Gross Weight*12 sets( Including out	7.6kg ( ± 5% )

packing )	
Quality guarantee	Main unit – 1 year
Manufacturing Date	See the packing

**Packing List:**

No.	Item	Qty.	Mark
1	Main unit	1	ProJet 2
2	Supragingival sandblasting handpiece	1	
3	Certificate	1	
4	Periodontal nozzle	1 pack	20 pcs in 1 pack
5	Long cleaning needle	1	
6	Short cleaning needle	2	
7	Instruction manual	1	
8	Powder chamber cap	1	
9	Sanblasting nozzle wrench	1	
10	Seal rings	1 Set	8 pcs in 1 pack
11	Periodontal sandblasting handpiece	1	
12	Cleaning brush	1	
13	Quick coupling	1	
14	Wrench	1	

**Specific please in kind prevail!**