




TAVOLA DI COMPARAZIONE

	VRN-Q6	DQ-80	DQ-40
Appearance			
Launch sales	2018	2023	2023
Control	3.5" LCD display	7" LCD touch screen	Touch panel
Handpiece	PH-1 LED Scaling Handpiece×2 AH-1 Supragingival Sandblasting Handpiece ×2 AH-2 Subgingival Sandblasting Handpiece AH-2 ×1	PH-1 LED Scaling Handpiece ×2 AH-1 Supragingival Sandblasting Handpiece ×2 AH-2 Subgingival Sandblasting Handpiece ×1	PH-1 LED Scaling Handpiece×1 HAH-1 Supragingival Sandblasting Handpiece ×1
Water Tank	700 ml ×2	900 ml ×1	900 ml ×1
Powder containers	2 pcs	2 pcs	1 pcs
Power	100~240V, 50/60Hz	110 or 230 V, 35 VA	110 or 230 V, 35 VA
Foot pedal	Multifunctional wireless foot pedal (IPX4 waterproof rating) ×1, Wired foot pedal (IPX1 waterproof rating) ×1	Multifunctional wireless foot pedal (IPX4 waterproof rating) ×1, Wired foot pedal (IPX1 waterproof rating) ×1	Multifunctional wireless foot pedal (IPX4 waterproof rating) ×1, Wired foot pedal (IPX1 waterproof rating) ×1
Tips	G1 ×2, G2 ×2, G4 ×2, G5 ×1, P1 ×2, P4 ×1, P11 ×2, P12 ×2, P12L ×2, P12R ×2, P16 ×2, IM1 ×1, E1 ×1 (Total 22 pcs)	G1 ×1, P1 ×1, G5 ×1, P11 ×1, P12 ×1, P12L ×1, P12R ×1, P16 ×1, IM1 ×1, E1 ×1 (Total 10 pcs)	G1 ×1, P1 ×1, G5 ×1, P11 ×1, P12 ×1, P12L ×1, P12R ×1, P16 ×1, IM1 ×1, E1 ×1 (Total 10 pcs)
Smart features	① Adjusting the amount of powder supplied at constant pressure ② Filter and check valve protect the powder blasting tip from clogging ③ Installed tips recognition function ④ Long press cleaning function	① Water heating system ② Adjusting the amount of powder supplied at constant pressure ③ Filter and check valve protect the powder blasting tip from clogging ④ Installed tip recognition function ⑤ Long press cleaning function	① Adjusting the amount of powder supplied at constant pressure ② Filter and check valve protect the powder blasting tip from clogging ③ Installed tip recognition function ④ Long press cleaning function
Other accessories	Metal torque wrench ×2 File holder key ×1 Box for sterilization ×2 + Tip care kit	Metal torque wrench ×2 File holder key ×1 Box for sterilization ×2 + Tip care kit	Metal torque wrench×1 File holder key ×1 Box for sterilization ×1 + Tip care kit