

R1 Plus

Gutta-percha Cutter



Designed by international dental experts

Conducting 220 precise manufacturing processes on the equipment imported from Germany and Japan




01

Gutta-percha point trimming,
smooth and glabrous cut surface



02

18 trimming and diameter measuring
holes with an accuracy of 0.01mm,
suitable for all gutta-percha points
on the market



03

Anti-rust and double-sided cutting
blade, smooth cutting can be done
by sliding forward or backward



04

30mm length measuring ruler



05

Laser marking guarantees durable
and fadeless marks



06


Made of aluminum alloy,
autoclavable, light weight, no rust



Guilin Woodpecker Medical Instrument Co., Ltd.
Information Industrial Park, National High-tech
Zone, Guilin, Guangxi, 541004 P. R. China.

Registrant/manufacturer/after-sales service provider:
Guilin Woodpecker Medical Instrument Co., Ltd.
Sales Dept.: 0773-5873196/2350599
Fax: 0773-5822450
After-sales Service Dept: 0773-5827898/13978361362
E-mail: woodpecker4@glwoodpecker.com
Website: <http://www.glwoodpecker.com>

Dealers



ZMN-XC-035